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


BYOVENT
THE FIRST INTENSIVE CARE
MECHANICAL VENTILATOR
OF TURKEY

biosys
MILLI MEDICAL SYSTEMLER

BYOVENT



1987

100%

Made In Türkiye

Biyovent is produced with the support of the Republic of Turkey Ministry of Industry and Technology. The Scientific and Technological Research Council of Turkey and Bilkent Cyberpark after a five-year research and development process.

The five-year design is verified and manufactured with the cooperation of engineers of defense industry and doctors who are experts at respiratory physiology.

BIYOVENT THE FIRST HIGH LEVEL INTENSIVE CARE MECHANICAL VENTILATOR OF TURKEY

DESIGNED AND PRODUCED FOR INTENSIVE CARE AND REANIMATION UNITS.

- Modern and ergonomic lines.
- User-friendly interface.
- Perfect performance.
- Traditional and innovative operating modes.
- Compatible with pediatric and adult patients.
- Low cost and maintenance.

Biyovent Respiratory Therapy Equipment


Ergonomic Design

- Design with smooth and modern lines
- Right/left and up/down angled, 15 inches, high resolution full touch screen monitor
- Touch-operated rapid access keys
- Top section detachable from support legs
- User-friendly expiration valve
- Protective carrying handles
- Shock-absorbing and impact resistant wheels
- External humidification support

Smart Safety System and User-friendly Interface


- Smart alarm identification and alarm silence (2min)
- Gradual auditory and visual alarm
- Adjustable apnea time and apnea backup mode (5-60 sec.)
- Automatic bilateral apnea ventilation mode
- 2 minute supply of oxygen (O2 suction)
- Stand-by mode
- Leak and trigger compensation
- Automatic tube compensation
- Comparative measurement of sensors and automatic calibration (when turning on the device and on request)

- Oxygen sensor
- Monitoring the trend of a patient for 1 week
- Logging the system for 6 weeks
- 2 hour internal battery
- 8 hour optional battery
- 5 ms valve response time
- 100 mbar emergency valve
- 50 mbar automatic expiration evacuation
- IP 21 impermeability
- Low air and oxygen pressure detection
- Automatic change of source
- Software update
- Working with a medical compressor
- Display of loss of main power and level of battery
- Display of technical failure, fan failure and connection loss alarm


Economic Solutions

- 100 % made in Türkiye
- Minimum maintenance cost
- Reusable expiration valve
- 2 years warranty + 3 years optional warranty extension
- Quick and qualitative technical service

Pediatric
and Adult

One Machine Compatible with All Patients

Biyovent can support all female & male and pediatric & adult patients.

It can deliver up to 150 breaths and go down to 20 cc low tidal.

It can calculate the values of patient elastance and compliance fast and precisely.

It can automatically calibrates itself when it is turned on.


Automatic tube compensation.

It provides leak compensation up to 80 %.

Inspiration hold and expiration hold in 1-60 sec. intervals.

Advanced
Adaptive Control

Biyovent uses advanced adaptive control algorithms. It responds in milliseconds.


Biyovent has Nasal Cpap and High Flow Oxygen Therapy modes.

Biyovent has an integrated nebulizer system.


Advanced
Adaptive Control

Biyovent works both in traditional and innovative modes.


- | | | |
|---------------------|-------------------|-----------------------|
| Pressure Controlled | Volume Controlled | Spontaneous and Smart |
| Modes: | Modes: | Modes: |
| P-ACV | V-ACV | SPN-PS |
| P-SIMV+PS | V-ACV(PRVC) | SPN-VS |
| P-CMV | V-CMV | |
| P-PSV | V-SIMV+PS | |
| P-Bilevel | V-SIMV(PRVC)+PS | |
| APRV | | |

Interface of Biyovent and Features of Software

Just identify the patient
in Biyovent for the rest.


Biyovent has finger compensation feature.

Technical Features of Biyovent

Types of Patients:

Pediatric and Adult

Modes of Ventilation

- P-ACV Pressure Controlled, Assisted Ventilation
- P-SIMV+PS Pressure Controlled, Synchronized Mandatory Ventilation With Pressure Support
- P-PSV Pressure Controlled, Ventilation With Pressure Support
- P-BILEVEL Pressure Controlled, Two-Level Ventilation
- P-CMV Pressure Controlled, Continuous Mandatory Ventilation
- APRV Airway Pressure Release Ventilation
- V-ACV Volume Controlled, Assisted Ventilation
- V-ACV(PRVC) Volume Targeted, Pressure Controlled, Assisted Ventilation
- V-CMV Volume controlled, Continuous Mandatory Ventilation
- V-SIMV+PS Volume controlled, Synchronized Mandatory Ventilation With Pressure Support
- V-SIMV(PRVC)+PS Volume Targeted, Pressure Controlled, Synchronized Mandatory Ventilation With Pressure Support


- SPN-PS Spontaneous Ventilation With Pressure Support
- SPN-VS Spontaneous Ventilation With Volume Support

nCPAP Nasal CPAP Mode

High Flow Oxygen Therapy Mode, 2-120 ltr/min
Spontaneous Breath Indicator

Compatibility with Patient and Performance

	Pediatric	Adult
Respiration rate	10-60 per minute	10-60 per minute
Inspiration	0-10 Sec	0-10 Sec
Tidal Volume	0-60 ml	0-120 ml
Flow	0-120 ltr per minute	0-120 ltr per minute
Trigger Sensitivity	0-120 ml	0-20 ml
Automatic Tube Compensation	0-80%	0-80%

Working Features

- Inspiration pressure 2-100 mBar
- Inspiration Time 0-1-10 sec
- Peep Pressure 1-50 mBar
- Respiratory rate (p): 1-150/min
(a): 1-100/min
- Tidal Volume (p): 20-600 mL
(a): 100-3000 mL
- Flow Rate (p): 1-60 ltr/min
(a): 1-120 ltr/min
- C2 Mixture 21-100%
- Spontaneous Pressure Support 0-100 mBar
I/E rate 1:10(x60*):10:1

(p): pediatric (a): adult

Detailed Features

- Apnea Time 5-60 sec, Bilateral Apnea Ventilation
- Apnea Mode P-CMV, V-CMV
- Flow Trigger 0.1-20 ltr/min
- Pressure Trigger Termination 0.1-20 mBar
- Termination Of Inspiration 0-80%
- Tramp 0.1-5 sec
- Automatic Tube Compensation 0-80%
- Automatic Leak Compensation 0-80%
- Trigger Compensation On-Off
- Inspiration Pause 1-60 sec
- Expiration Pause 1-60 sec
- Internal (Integrated) Nebulizer 1-20 Ltr/min
- O2 Support 2 min

Displayed Data on screen

P Peak	Measurement of Peak Pressure Inspiration	WOB	Energy Spent During Inspiration
P Peep	Measurement of PEEP	WOB/Lt	Energy Spent During Inspiration / Volume
P Plateau	Measurement of Plateau Pressure Inspiration	V Residual	Residual volume at the end of breathing
P Average	Measurement of Average pressure	V Expiration	Expiration Tidal Volume
F Inspiration	Inspiration Flow	V Ads	Anatomic Dead Space Measurement
F Expiration	Expiration Flow	AutoPeep	Trapped Air Pressure After Respiration Occlusion
MVE	Volume Measurement	PO1	Pressure Measurement Per 100 ml/Seconds
SPM/MVE	Spontaneous Minute Volume Measurement	RSBI	Rapid Shallow Breathing Index
SPM/Ve/MVE	Spontaneous Volume Per Minute / Volume Ratio Per Minute	PTP	Negative Pressure x Negative Pressure Time
V Tidal	Tidal Volume	FTP	Negative Flow x Negative Flow Time
FI02	Oxygen Ratio	P NIF	Negative Inspiration Pressure Force
Respiratory rate	Number of Breaths Per Minute	MV/spk	Spontaneous Ratio to Mandatory Minute Volume
Spontaneous respiratory rate	Spontaneous Number of Breaths Per Minute	Leak Rate	Leak Volume Rate
T Inspiration	Inspiration Time	Leak Volume	Leak Volume After Respiration Cycle
T Expiration	Expiration Time		
I/E	Inspiration and Expiration Time Ratio		
R Airway	Airway Resistance Measurement		
C Static	Static Compliance Measurement		
C Dynamic	Dynamic Compliance Measurement		
Elastance	Elastance Measurement		
RC Constant	RC Time Constant Measurement		

Alarm Features

Auditory and Visual Alarm and Recording
Two-Minute Alarm Silence

Inspiration Pressure Lower Limit / Upper Limit
Tidal Volume Lower Limit / Upper Limit
Speed of Respiration Lower Limit / Upper Limit
Volume Per Minute Lower Limit / Upper Limit
I/E Ratio Lower Limit / Upper Limit
FI02 Lower Limit / Upper Limit
Apnea Time Upper Limit
Leak Upper Limit

Graphic Features

Pressure Time Graphic
Flow Time Graphic
Volume Time Graphic
Pressure-Volume, Pressure-Flow, Volume-Flow Cycles

Patient Records and Logging

Last 3 Days Record of Ventilation Value Trend and Graphical Representation.
Information Storage With The Capacity of 2000 Record Alarm and Ventilation Alarm.

Electrical Features

Battery Time 2 Hours + 8 Hours Optional
Mains Voltage Power 180 - 264 VAC
Consumption 47-63 Hz 100W

Features of The Source of Pressure

O2 Pressure 2.5 - 7 Bar Central System / Tube
Air Pressure 2.5 - 7 Bar Central System / Tube
Automatic Change and Alarm Display When The Source Is Consumed
Working With a Medical Compressor or Regulator

Size and Weight

Length 150cm
Depth 44cm
Width 42cm
Weight 55kg
Monitor Movement Left and Right 150°
Up and Down 15°
15 inch Full Touch Screen Monitor
Pendant and Column Mounting

Digital Interfaces

4 USB, 2 COM, 2 Ethernet

Comparison of Modes

Biosys	Puritan Bennett	Dräger	GE	Hamilton	Maquet	Mindray
P-ACV	A/C: PC	PC-AC	PCV	PCV	PC	P-AC
P-SIMV+PS	SIMV: PC	PC-SIMV	SIMV-PC	PSIMV+	SIMV-PC+PS	P-SIMV
P-PSV	PS	SPN-CPAP/PS	CPAP/PSV	Spont	PS	PSV
P-BiLevel	BiLevel	PC-BIPAP	BiLevel	DuoPAP	Bi Vent	DuoLevel
APRV	APRV	PC-APRV	APRV	APRV	BiVent-APRV	APRV
V-ACV	A/C: VC	VC-AC	VCV	(S)CMV	VC	V-AC
V-ACV(PRVC)	VC	AutoFlow	PCV-VG	APV/SIMV	PRVC	PRVC
V-CMV	A/C: VC	VC-CMV	VCV	CMV	VC	V-AC
V-SIMV+PS	SIMV-VC	VC-SIMV	SIMV-VC	SIMV	SIMV-VC+PS	V-SIMV
V-SIMV(PRVC)+PS	VC+	VC-SIMV+ AutoFlow	SIMV-PCVG	APV/SIMV+	SIMV-PRVC+PS	PRVC
SPN-PS	PS	SPN-CPAP/PS	CPAP	Spont	PS/CPAP	-
SPN-VS	VS	SPN-CPAP/VS	-	-	VS	-

- +PS(Pressure Support) feature supports the breathing efforts with pressure.
- PRVC(Pressure Regulated Volume Control) feature provides pressure control for volume target.

CERTIFICATE

Full Quality Assurance System

Medical Devices Directive 93/42/EEC Annex II (Excluding Section 4)

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Related Directives and Annex : 93/42/EEC Medical Devices Directive - Annex II (Excluding Section 4)

Product : Non-Sterile-Intensive Care Ventilator - Class IIB

GMDN : 47244

Certificate Number : M.2020.106.13453
Report Number : MD.4068.IB
Initial Assessment Date : 14.02.2020
Registration Date : 30.03.2020
Revision Date /No : -
Expiry Date : 27.05.2024

UDEM International Certification
Auditing Training Centre Industry
and Trade Inc. Co.

UDEM hereby declares that the requirements of Annex II, excluding section 4 of the 93/42/EEC Directive have been met for the listed products. The above named manufacturer has established and applied a quality assurance system, which is subject to periodic surveillance audits, defined by Annex II, section 5 of the listed standard directive. According to Annex II, section 4 an EC design examination certificate has been issued for the Class II devices on the market. UDEM responsibility for class II devices covered by the EC certificate is limited to manufacturing have related to its reporting and not related to field conditions. The above stated manufacturer has submitted to product conformity with technological requirements. UDEM measurement laboratory, its certificate number at the proposal of UDEM, Verification of Certification Auditing Centre, Ankara, Turkey, has issued a certificate for the product. UDEM hereby declares that the manufacturer is responsible for the validity of the certificate for 5 years from the registration of the certificate. UDEM will not renew the validity of the certificate if the manufacturer with the approved product to UDEM. UDEM will not renew the validity of the certificate in question, the mentioned company should stop placing the product on the market. The validity of the certificate can be checked through www.udem.com.tr.

Address: Mülkent Mahallesi 2073 Sokak, (Eski 93 Sokak) No:10 Çankaya - Ankara - TURKEY
Phone: +90 0312 443 03 90 Fax: +90 0312 443 03 74
E-mail: info@udem.com.tr www.udem.com.tr

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LVT Test Laboratuvarları Ltd. Şti. accredited by TÜRKAK under registration number AB-0341-T for IEC/ISO/TS/EN 17025:2017 as test laboratory.

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The Turkish Accreditation Agency (TÜRKAK) is a signatory to the European Cooperation for Accreditation (EA) Multilateral Agreement with the International Laboratory Cooperation (ILCO) Mutual Recognition Arrangement (MRA) for the recognition of test reports.

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Mühür

Deney Sorumlusu
Person in Charge of Test

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
This report shall not be reproduced either in full or in part without the written permission of the laboratory.
This report shall not be reproduced either in full or in part without the written permission of the laboratory.

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Beyan Frekans : 50 Hz

Beyan Koruma Derecesi : 20

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Çihaz - Malzeme Listesi : Class I and Internally powered ME equipment.

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Deneysel Bilgisi Test details	Çihaz kodu Device code	Ölçülen değer Measured value	Ölçüm belirsizliği Measurement uncertainty
Ambient Temperature	LC349	Cl. 3.1	%1.6
Ambient Moisture	LC349	Cl. 3.2	±2.98
Creepage distances and air clearances	LC365	Ins. Diagram	±0.0758
Creepage distances and air clearances	LC443	Ins. Diagram	±0.06
Power input	LC96	See Table 4.11	±1.31
Humidity preconditioning treatment (Temperature)	LC215	See Clause 5.7	±1.26
Humidity preconditioning treatment (Humidity)	LC215	See Clause 5.7	±3.91
Accessible parts and applied parts	LC30	See Table 8.4.2	±1.34
ME equipment intended to be connected to a power source by a plug	LC30	See Table 8.4.3	±1.34
Impedance and current-carrying capability	LC85	See Table 8.6.4	±3.96
Leakage currents and patient auxiliary currents	LC91	See Table 8.7	±3.6
Dielectric strength	LC85	See Table 8.8.3	±1.91
Mechanical strength and resistance to heat	LC100	See Table 8.3.4.1	±2.67
Gaps Measurement	LC365	See Table 9.2.2.2	±0.0539
Gaps Measurement	LC443	See Table 9.2.2.2	±0.06
Instability in transport position	LC284	See Table 9.4.2.1	±1.33
Instability excluding transport position	LC284	See Table 9.4.2.2	±0.66
Audible acoustic energy	LC44	See Clause 9.6.2.1	±0.31
ME equipment not intended to produce diagnostic or therapeutic x-radiation	LC320	See Table 10.1.1	±0.0612
Maximum temperature during normal use	LC74	See Table 11.1.1	±1.39
Push test	LC204	See Table 15.3	±0.59
Mechanical stress relief test	LC100	See Table 15.3	±2.67

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Numune Sample	Uygulanan Deneysel Test Applied Test	Sonuç Result
Biyovent	TS:EN 60601-1:2009+A1:2014+AC:2011 #AYZ:2015+A1/AC:2014 IEC 60601-1:2005/AMDT:2012/CER1:2014	OLUMLU Passant

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3 (75±3) %RH

Deneysel Metodundan

Sapma, Ekleme ve Çıkarımlar
Deviations, Additions & Subtractions from the Test Method

Deneysel standart deneysel metoduna göre uygulanmıştır.
Tests were made according to the clauses of the relevant standards

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Medical Electrical Equipment Tests

9. Deneysel Uygulamalar:

Test Applications

IEC 60601-1

Medical electrical equipment

Part 1: General requirements for basic safety and essential performance

Test specification:

Standard IEC 60601-1:2005, COR1:2006, COR2:2007, AMB1:2012
(or IEC 60601-1:2012 reprint)

Test procedure Type Test

Non-standard test method.....

Test Report Form No IEC60601_1P

Test Report Form Originator UL(US)

Master TRF 2019-10-11

Copyright © 2019 Worldwide System for Conformity Testing and Certification of Electrotechnical Equipment and Components (IECEE), Geneva, Switzerland. All rights reserved.

This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as the copyright owner and notice of the material, IECEE, takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material for its placement and context.

If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the references to the CB Scheme procedure shall be removed.

This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IEC60601.

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Medical Electrical Equipment Tests

©copy or marking plates

The artwork below may be only a draft. The use of certification marks on a product must be authorized by the respective NCBs that own these marks.

Biopsy Biomedical Malhenderlik-San. Tic. Ltd. Sti.
ATG Bolgesi, Gopnarak, Cihangirpazari A Blok, 303C
Cankaya/Beşiktaş/İstanbul, Turkey

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Medical Electrical Equipment Tests

GENERAL INFORMATION

Test item particulars (see also Clause 6):	Mobile
Classification of installation and use:	Equipment
Device type (component/sub-assembly/equipment/system):	Intensiv care unit
Intended use (including type of patient, application location):	Continuities
Mode of operation:	Internally powered / appliance-coupler
Supply connection:	Control separate/dedicated.
Accessories and detachable parts included:	None
Other options included:	None
Name and address of factory(ies):	ARGELIK A.Ş. Elektronik İşletmesi: Çarşakçy Organize Sanayi Bölgesi Karagağ Mah. 3. Sokak No:1A 59510 Kapaklı/Tekirdağ
Testing	
Date of receipt of test item(s):	16.04.2020
Dates tests performed:	16.04.2020 - 28.04.2020
Possible test case verdicts:	N/A
- test case does not apply to the test object:	Pass (P)
- test object does meet the requirement:	N/E (collateral standards only)
- test object was not evaluated for the requirement:	Fail (F)
- test object does not meet the requirement:	
Abbreviations used in the report:	
- normal condition:	N.C.
- means of Operator protection:	M.OOP
- single fault condition:	S.F.C.
- means of Patient protection:	M.OBP
General remarks:	

Before starting to use the TRF please read carefully the 4 instructions pages at the end of the report or how to complete the new version "K" of TRF for IEC for 60801-1 3rd edition with Amendment 1.

"See Attachment #)" refers to additional information appended to the report.

"See appended table)" refers to a table appended to the report.

The tests results presented in this report relate only to the object tested.

This report shall not be reproduced except in full without the written approval of the testing laboratory.

List of test equipment must be kept on file and parallel to the test report.

Additional test data and/or information provided in the attachments for this report.

Throughout this report a comma / point is used as the decimal separator.

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